



ABLE LIFE™
Comfortably Independent

EZ DOOR KNOB GRIPS

OPEN ANY
DOOR WITH
THE TOUCH
OF YOUR
FINGER



OPENING DOORS IS EASY AGAIN

PRODUCT # 8130

EXPERIENCE ABLE LIFE



Open doors with the touch
of your finger



Anti-bacterial material



Glow-in-the-dark inserts

FEATURES & BENEFITS

- Low cost alternative for door handles
- Stretches over most doorknobs
- Anti-bacterial material
- Removable glow-in-the-dark inserts
- Keyhole easily accessible

SPECIFICATIONS

- **Door Types:** Fits most round door knobs
- **Grip Width:** 3"
- **Product Dimensions:** 3" x 2" - 1.6" diameter 0.2 lbs
- **Package Dimensions:** 7" x 2.75" x 2" 0.2 lbs
- **Case of 24 Dimensions:** 9" x 9" x 9" 5 lbs
- **Made Of:** Silicone

ADA Compliant



EZ Door Knob Grips Slip on Quickly- No Tools Needed!



Moisten the Grips with water



Stretch over door knob



Continue to work onto
door knob until covered